Name: Date:

## Equilibrium Potentials I Lab Notebook Reading: Equilibrium Potentials I

Purpose of Lab:
Problem Set 1
<b>a.</b> If you have not already done so, download the Nernst Potential Simulator using the links in Canvas.
<b>b.</b> Double click the program icon to open the simulator. Under the File menu, choose Load Initial Conditions, and then select <i>Problem1.init</i> . This file should be in the same folder as the one in which you found the Nernst Potential Simulator program. Note that the <b>permeability</b> , which is a measure of the number of channels that will allow a particular kind of ion to cross the membrane, has been set to zero for all of the ions. When the permeability is set to 1, the maximum number of channels are available to those ions. Also note that the <b>Electrostatics</b> box is unchecked. This means that the ions will be treated as uncharged particles. Consequently, <b>ignore the Membrane Potential plot</b> on the right side of the window. It is only applicable when Electrostatics is turned on, which will only happen during part (f) of this problem.
<b>c.</b> Start the simulation (press Start). Click on one of the red ions, which will enlarge it. <i>Describe</i> what you observe about the motion of the single ion, as well as the movements of all of the ions in the intracellular compartment. Predict the long term behavior of the system, record your prediction, and observe the system for at least 2,000 iterations to see if your long term prediction is likely to be correct (the iteration number is updated on the lower left-hand corner of the window).
d. Pause the simulation (press Pause). Change the permeability for the K+ ions from 0 to 1 (either by moving the slider, or typing 1 into the box). Also change the permeability of the Na+ and Cl- ions to 1. Note that the Selective Permeability box is unchecked. Thus, increasing the permeability to Na+ and Cl- under these conditions will provide many opportunities for potassium ions to move through the membrane. You can see the many channels shown in the zoom window at the bottom middle of the simulation. <b>Before continuing, take a screenshot of the simulation for your lab notebook.</b>

	Paste screenshot here. Make sure to label the figure, you can do this by right clicking the image and clicking "Insert Caption"	
you have sele and extracellul final intracellul in the simulation	the simulation (press Continue). Observe what happens to the individual ected, and to all of the ions over time. Also observe the change in the in allar concentrations of K+ listed in the table under the graph to the right. Follar and extracellular concentrations of K+ ions. Note that in this exercise ion have no charge (the Electrostatics setting is disabled), and so they welly-neutral particles, such as glucose. Run the system for up to another.	tracellular Predict the the "ions" rill behave

the simulation for your lab notebook.

Paste screenshot here. Make sure to label the figure, you can do this by right clicking the image and clicking "Insert Caption"

iterations to see if your long term prediction is likely to be correct. Take another screenshot of

	ompartment. Change the permeability of all the ions from 0 to 1. <b>Before starting 1, take another screenshot for your lab notebook.</b>
	Paste screenshot here. Make sure to label the figure, you can do this by right clicking the image and clicking "Insert Caption"
extracellular configuration if your long tention compartments identical volumes.	ation, and again click on one of the red ions. Predict the long term intracellular and oncentrations of K+ ions, and observe the system for about 5,000 iterations to see in prediction is correct. Note to the careful student: the extracellular and intracellular have slightly different volumes in this simulation, but you may treat them as nes for these questions. Take another screenshot of the simulation after at least ins have passed.
	Paste screenshot here. Make sure to label the figure, you can do this

by right clicking the image and clicking "Insert Caption"

**e.** Take the simulation back to the initial conditions loaded from the file by pressing the "Reset with Loaded Initial Conditions" button. Now switch the intracellular concentration of K+ from 2000 to 0, and the extracellular concentration of K+ from 0 to 2000, so that the K+ ions are all in the

Please make sure to do the reading prior to completing the problem set.				
species has different concentrations in the intra- to freely diffuse across a permeable membrane uncharged (as in the exercises above). Then, dis (e.g., is an ion) and is associated with a second and chloride ions, or potassium and chloride ion concentrations as the first on the same side of the prediction by running the simulation and see who	ehavior of any system in which a single particle cellular and extracellular spaces, and is allowed a. First answer this assuming that the particle is cuss what will happen if the particle has a charge in it is it is in of opposite and equal charge (e.g., sodium as); assume the second ion has the same initial emembrane and is equally permeable. Test your at happens. You can give the particles charge in the contract of the contr			
ste screenshot here. Make sure to label the figure, you can do this by right clicking the image and clicking "Insert Caption"	Paste screenshot here. Make sure to label the figure, you can do this by right clicking the image and clicking "Insert Caption"			

## Question 2: What about osmotic balance?

You must follow the instruction in the reading for this problem. By the first rule, water should move from the side that has more water molecules (the left side) to the side that has less (the right side), shrinking the size of the cell. Using **Equation 3**, and assuming for simplicity that each volume is 1 liter initially, it is easy to show that that the volume change is 9/11 liter. *Please do the calculation to make sure you see that this is true! Please include your work and the results of your calculation in your notebook.* Using this result in **Equation 2**, it is easy to show that the total concentration of particles on each side of the membrane will now be equal.

**Question 3:** What happens to the Nernst potential for chloride under those circumstances?

You must follow the instruction in the reading for this problem. Examine Equation 5. The concentration of the chloride ion will now be the same inside and out. What will the resulting voltage be? *Please do this calculation. Please include your work and the results of your calculation in your notebook.* If you don't find it easy, please review your understanding of logarithms.

**Question 4:** Calculate the concentrations of ions and equilibrium across the membrane.




Figure 6: Balancing osmotic pressure, bulk electroneutrality, and concentration and electrical gradients

You must follow the instruction in the reading for this problem. We need to ensure bulk electroneutrality within the cell, that is, there have to be as many positive as negative charges inside. What is the internal concentration of chloride ions based on this rule (i.e., the concentration of chloride ions on the left side of Figure 6)?

Since we've filled in the concentrations of all components within the cell, what is the total concentration of solute particles inside the cell (left side of Figure 6)?
We have to maintain bulk electroneutrality outside of the cell as well; additionally, to guarantee osmotic balance, we have to ensure that the concentration of solute particles outside the cell equals the concentration of solute particles inside the cell. Given this information, what are the concentrations of the sodium and chloride ions outside the cell (right side of Figure 6)?
Now that you have the concentrations of chloride ions both outside and inside the cell, what is the equilibrium or Nernst potential across the cell membrane, i.e., please calculate $E_{cl}$ using Equation 5.
Once you've tried the calculations, you can check the results in the reading.
Discussion:
•
•